**Installation Instructions | Installationsanleitung | Manuel d'installation |
Manuale di installazione | Instrucciones de instalación | Instruções de instalação**
Translation of the Original Installation Instructions | Original-Installationsanleitung |
Traduction du manuel d'installation original | Traduzione del manuale originale |
Traducción de las instrucciones de instalación originales | Tradução das instruções de
funcionamento originais

MD8 Airscan Sampling Head | MD8 Airscan Sammelkopf | Tête de prélèvement MD8 Airscan | Gruppo aspirante MD8 Airscan | Unidad de aspiración MD8 Airscan | Unidade de aspiração MD8 Airscan

16746-01--SHB | 16746-01--SHT

Sampling Head for Airborne Microbe Collector | Sammelkopf für Luftkeimsammler |
Tête de prélèvement pour collecteur de germes | Gruppo aspirante per campionatori
d'aria | Unidad de aspiración para colector de gérmenes aéreos | Unidade de
aspiração para amostradores de ar

English – page 3

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Contents

1 About these Instructions	4	4 Installation	12
1.1 Validity.....	4	4.1 Equipment Supplied.....	12
1.2 Accompanying Documents	4	4.2 Prerequisites for Installation ..	12
1.3 Symbols Used	4	4.3 Unpacking the Equipment.....	13
1.3.1 Warnings.....	4	4.4 Setting up the Sampling Head	13
1.3.2 Other Symbols.....	4	4.5 Conducting Exhaust Air Away	17
2 Safety Instructions	5	5 Technical Data	19
2.1 Intended Use.....	5	5.1 General Data	19
2.2 Significance of These Instructions	6	5.2 Ambient Conditions	19
2.3 Working on the Electrical Equipment.....	6	6 Accessories and Consumables ...	20
2.4 Accessories, Consumables and Spare Parts.....	6	6.1 Accessories.....	20
3 Device Description	7	6.2 Consumables	21
3.1 Device Overview	7		
3.2 Bayonet Connection.....	8		
3.3 Tri-Clamp Connector (TC 50.5).....	8		
3.4 Air Inlet	9		
3.5 Air Outlet.....	10		
3.6 Adapter Bayonet on Tri-Clamp TC 50.5.....	10		
3.7 Filter Adapter	11		
3.8 Airborne Microbe Collection ..	11		

1 About these Instructions

1.1 Validity

These instructions apply to the device in the following versions:

Device	Type
Sampling head with bayonet connection	16746-01--SHB
Sampling head with Tri-Clamp connection (TC 50,5)	16746-01--SHT

1.2 Accompanying Documents

- ▶ In addition to these instructions, observe the following documents:
 - MD8 Airscan Command Unit Operating Instructions
 - Documentation on the filter
 - If the sampling head is installed in a system: Instructions for the system

1.3 Symbols Used

1.3.1 Warnings

NOTICE

Denotes a danger that can result in property damage if the risk is **not** avoided.

1.3.2 Other Symbols

- ▶ Required action: Describes actions that must be performed.
- ▷ Result: Describes the result of the actions carried out.

2 Safety Instructions

2.1 Intended Use

The MD8 Airscan Sampling Head serves as a vacuum source to pump ambient air. The sampling head is suitable for the following applications:

- Use in controlled areas such as clean rooms
- Use in systems, e.g. sterility isolators, filling systems, blow-fill-seal machines
- Use in non-controlled areas such as laboratories

The sampling head must be operated with the MD8 Airscan Command Unit and with a filter. The sampling head **cannot** be operated without a command unit. The filter is designed for single use and must be disposed of each time after operation.

Do **not** allow the sampling head to come into contact with liquids. Do **not** use the sampling head in potentially explosive atmospheres. Only use the sampling head indoors.

Only use the sampling head with the equipment and under the operating conditions described in the Technical Data. Do **not** rebuild or technically alter the sampling head.

Any further use beyond this is considered **improper** and is prohibited. These instructions are part of the product. The device is exclusively designed for use according to these instructions.

2.2 Significance of These Instructions

Failure to follow the instructions in this manual can have serious consequences, e.g. exposure of individuals to electrical, mechanical or chemical hazards.

- ▶ Before working with the device: Read the instructions carefully and completely.
- ▶ If the instructions are lost: Request a replacement or download the latest version from the Sartorius website (www.sartorius.com).
- ▶ The information contained in these instructions must be made available to all individuals working with the device.

2.3 Working on the Electrical Equipment

Any work on or modifications to the electrical equipment of the device may only be carried out by Sartorius Service personnel. The device may only be opened by Sartorius Service personnel.

2.4 Accessories, Consumables and Spare Parts

The use of unsuitable accessories, consumables and spare parts can affect the functionality and safety of the device and have the following consequences:

- Risk of injury to persons
 - Damage to the device
 - Malfunctions of the device
 - Device failure
-
- ▶ Only use accessories, consumables and spare parts from Sartorius (for accessories, see the MD8 Airscan Command Unit Operating Instructions). Information on operational quality is available upon request from Sartorius.
 - ▶ Only use accessories, consumables and spare parts that are in technically perfect condition.

3 Device Description

3.1 Device Overview


Fig. 1: Sampling head with bayonet connection (example)

No.	Designation	Description
1	Filter holder	Is designed as a bayonet connection or as a Tri-Clamp connection (TC 50.5).
2	Air outlet	Is designed as a Tri-Clamp connector (TC 50.5).
3	Connection for command unit	Accommodates the command unit connection cable.
4	Device base (ring-shaped)	Protects the sampling head against slipping. Not depicted.
5	Manufacturer's ID label	Is affixed to the base plate. Not depicted.
6	Sampling head	Includes the device motor, fan and control-relevant components.
7	Air Inlet	

3.2 Bayonet Connection

The bayonet connector accommodates the filter adapter. Sampling heads with a bayonet connection are particularly suitable for applications where the sampling head is **not** installed in a system.

Fig.2: Bayonet Connection

No.	Designation	Description
1	Mounting	Accommodates the filter adapter.
2	Groove	Accommodates the seal. Inserted at factory.

3.3 Tri-Clamp Connector (TC 50.5)

Sampling heads with a Tri-Clamp connection (TC 50.5) are particularly suitable for installation in systems, e.g. isolators, filling systems and blow-fill-seal machines.

Fig.3: Tri-Clamp connector (TC 50.5)

No.	Designation	Description
1	Mounting	Allows installation in pipelines using a Tri-Clamp connection.
2	Groove	Accommodates the seal for the Tri-Clamp connection.

3.4 Air Inlet

The air to be filtered is drawn in at the air inlet.

Fig. 4: Air inlet at the sampling head with bayonet connection (example)

No.	Designation	Description
1	Stainless steel strainer	Prevents larger objects from entering the sampling head. Is installed inside the sampling head.
2	Sensor	Detects the amount and speed of air drawn in. Is installed inside the sampling head. Not depicted.
3	Gasket	Is installed inside the sampling head. Not depicted.

3.5 Air Outlet

At the air outlet, the air from the sampling head is directed into the ambient air or into pipes. The air outlet is designed as a connector piece. The connector piece can accommodate pipes and adapters.


Fig.5: Air outlet (example)

No.	Designation	Description
1	Connector piece	Allows installation in pipelines using a Tri-Clamp connection. Is designed as a Tri-Clamp 50.5.
2	Groove	Accommodates the seal for the Tri-Clamp connection.
3	Gasket	Is installed inside the connector piece. Not depicted.

3.6 Adapter Bayonet on Tri-Clamp TC 50.5

Sampling heads with a bayonet connection are particularly suitable for applications where the sampling head is **not** meant to be installed in a system. This allows for flexible operating conditions.

Sampling heads with a Tri-Clamp connection (TC 50.5) are particularly suitable for installation in systems, e.g. isolators, filling systems and blow-fill-seal machines. The Tri-Clamp connection (TC 50.5) allows the sampling heads to be connected to the pipes of the existing systems.

The adapter from Bayonet to Tri-Clamp TC 50.5 converts the bayonet connection into a Tri-Clamp connection. As a result, sampling heads with a bayonet connection can also be connected to the pipelines of systems.

3.7 Filter Adapter

The filter adapter accommodates the filter and, depending on the model, is suitable for accommodating gelatine membrane filters or cellulose nitrate membrane filters.

Fig. 6: Filter adapter for gelatine membrane filters

No.	Designation	Description
1	Filter mount	Accommodates the gelatine membrane filter.
2	Attachment for bayonet	Is placed on the bayonet connection.

The filter adapter must be placed on a bayonet connection. An adapter with bayonet connection is additionally required for sampling heads with a Tri-Clamp connection (TC 50.5).

3.8 Airborne Microbe Collection

Pumping of Air Sample Volume (Isokinetic Sampling)

The sampling head generates a vacuum and pumps a defined air sample volume from the ambient air. The air flow can be adapted to the ambient speeds so that the flow velocities in critical areas, e.g. in clean rooms, are **not** disturbed and turbulences are avoided.

Retention of Biological Particles

Active airborne microbe collection takes place in the filter. The biological particles from the air, e.g. viruses, bacteria, yeasts and fungi, are retained in the filter. This retention is almost independent of the configured air flow.

Incubation of the Filter

Since the filter does **not** contain germ-growth-influencing substances, the filter must be transferred to a suitable nutrient medium and incubated, e.g. to an agar plate with Sabouraud, for the detection of fungi and TSA for total microbial count determination.

4 Installation

4.1 Equipment Supplied

Product	Quantity
Sampling head with bayonet connection or Tri-Clamp connection (TC 50.5)	1
Calibration report	1
Installation Instructions	1

4.2 Prerequisites for Installation

Procedure

► Make sure that the following conditions are met at the place of installation:

Condition	Features
Ambient conditions	– Suitability tested (for ambient conditions, see Chapter 5.2, page 19)
Setup surface	– Stable and level surface – Sufficient room for the sampling head (for information on space requirements, see Chapter "5.1 General Data," page 19) – Sufficient bearing capacity for the sampling head (for weight, see Chapter "5.1 General Data," page 19).

4.3 Unpacking the Equipment

Procedure

- ▶ Unpack the device.
- ▶ If the device is stored temporarily: Observe the storage instructions (for storage information, refer to the MD8 Airscan Command Unit Operating Instructions).

4.4 Setting up the Sampling Head

Installation in controlled and non-controlled areas

The sampling head can be set up in non-controlled and controlled areas (controlled areas, see Chapter "5.2 Ambient Conditions," page 19).

When the sampling head is set up in controlled areas: The flow speeds in the controlled area can be affected by the operation of the sampling head and turbulences can occur. The sampling head suction rate must be adjusted to the flow conditions in the controlled area in order to allow for isokinetic sampling (to adjust the air flow, see MD8 Airscan Command Unit Operating Instructions).

In filling systems, isolators and blow-fill-seal systems, it is recommended that the sampling head be installed outside this critical area and that the filter be connected to the sampling head air inlet by means of pipes or hoses.

Sampling Head in Non-controlled Areas


Fig. 7: Sampling head set up in the non-controlled area, connected by a PVC hose (example)

No.	Designation	Description
1	Adapter or connector for Tri-Clamp connection	Connects the Tri-Clamp connection of the sampling head to silicone hose.
2	Hose	PVC hose
3	Hose adapter	PVC hose adapter
4	Sterile barrier	For example: Wall of a sterile isolator
5	Sterile filters	For example: Sartofluor® capsule on both sides with TC 50.5
6	Angle connector	Connects Tri-Clamp (TC 50.5) adapter to bayonet and sterile filter.
7	Tri-Clamp (TC 50.5) adapter on bayonet	Mount with bayonet connection
8	Filter adapter	Type: Adapter for gelatine membrane filter (GMF)
9	Filter	Type: Gelatine membrane filter (GMF)

Sampling Head in Controlled Areas


Fig.8: Sampling head set up in a controlled area, connected by a PVC hose; exhaust channeled into non-controlled area (example)

Pos.	Designation	Description
1	Exhaust line	PVC or silicone hose
2	Adapter	TC 50.5 adapter on hose nozzle
3	Sterile barrier	For example: Wall of a sterile isolator
4	TC 50.5 adapter on single-level hose nozzle DN 25 adapter	Is inserted in the air outlet.
5	Adapter or connector for Tri-Clamp connection	Connects the Tri-Clamp connection of the sampling head to silicone hose.
6	Hose	PVC hose
7	Hose adapter	PVC hose adapter
8	Sterile filters	For example: Sartofluor® capsule on both sides with TC 50.5
9	Angle connector	Connects Tri-Clamp (TC 50.5) adapter to bayonet and sterile filter.
10	Tri-Clamp (TC 50.5) adapter on bayonet	Mount with bayonet connection
11	Filter adapter	Type: Adapter for gelatine membrane filters (GMF)
12	Filter	Type: Gelatine membrane filter (GMF)

Setting up the Sampling Head in a Controlled or Non-controlled Area

Procedure

- ▶ Set up the sampling head at the intended installation location in the controlled or non-controlled area.
- ▶ When using a sampling head with a TC 50.5 Tri-Clamp connection: Mount a Tri-Clamp bayonet adapter on the sampling head (1) (see Chapter "6.1 Accessories," page 20).
- ▶ Place a filter adapter on the bayonet connector (1).



Installing Filter Adapter in Controlled Area

Requirements

- The sampling head is installed in a non-controlled area.
- The filter adapter is **not** mounted.

Procedure

- ▶ Install the filter adapter (1) in the controlled area. Use suitable accessories and consumables (2), e.g. connectors and sterile filters (see Chapter "6 Accessories and Consumables," page 20).
- ▶ Connect a hose or a pipe (2). Use suitable accessories (1), e.g. connecting pieces (see Chapter "6 Accessories and Consumables," page 20).

- ▶ If a sampling head with a bayonet connector is used: Mount bayonet to Tri-Clamp (TC 50.5) adapter on the bayonet connection.
- ▶ Insert a gasket for the Tri-Clamp connection in the groove (1) of the connector piece.
- ▶ **NOTICE** Damage to the sensor when connecting pipes or adapters! Insert the piping and adapters only so far as to **not** allow the stainless steel strainer and the sensor to be touched.
- ▶ Connect the hose or piping to the air inlet of the sampling head. Use suitable accessories, e.g. connectors, adapters or Tri-Clamps (see Chapter "6.1 Accessories," page 20).
- ▶ Secure the connection using a Tri-Clamp.

4.5 Conducting Exhaust Air Away


The exhaust air that flows out of the sampling head air outlet can be left in the ambient air or conducted into an exhaust air system.

If the sampling head is operated in controlled areas of clean room class A-B: The exhaust air should be removed.

When the sampling head is operated in controlled areas of clean room class C-D: It is recommended that the exhaust air be removed if the particle concentration in the exhaust air exceeds the permissible air purity limit. To this end, observe the applicable regulations on air purity, e.g. the valid United States Pharmacopeia (USP), Eu.Ph or ISO guidelines.

Procedure

- ▶ When the exhaust air is left in the ambient air:
Check the required clearance for the air outlet (for information on clearance, see Chapter "5 Technical Data," page 19).
 - ▶ If required: Position the sampling head differently.
 - ▷ The exhaust air flows into the ambient air.
- ▶ If the exhaust air is to be conducted away:
 - ▶ Insert a gasket for the Tri-Clamp connection in the groove (1) of the connector piece.
 - ▶ Connect a hose or a pipe to the connector piece (2). Use suitable accessories, e.g. connectors, adapters or Tri-Clamps (see Chapter "6.1 Accessories," page 20).
 - ▶ Secure the connection using a Tri-Clamp.

5 Technical Data

5.1 General Data

	Unit	Sampling head with bayonet connection	Sampling head with Tri-Clamp connection (TC 50.5)
Dimensions (L x W x H)	mm	128 x 112 x 137	128 x 112 x 154
Weight	kg	2.6	2.6
Power supply			
Only using MD8 Airscan Command Unit			
Material			
Seals		Silicone	Silicone
Surfaces		Stainless steel	Stainless steel

5.2 Ambient Conditions

	Unit	Value
Installation location		
Standard laboratory rooms and clean rooms according to ISO 5 to ISO 9;	m above sea level	< 3000
Use in controlled areas: Clean room class A-D		
Air outlet required clearance	mm	> 100
Recommended distance between sampling head and persons at the place of installation of the sampling head	mm	> 1000

	Unit	Value
Ambient temperature		
Operation	°C	4 – 40
Transport and storage	°C	4 – 45
Temperature of the motor in the sampling head, max.	°C	65
Relative humidity, max.		
Operation	%	80
Transport and storage	%	90
Electromagnetic compatibility according to DIN EN 61326-1: Electrical equipment for measurement, control and laboratory use – EMC requirements – Part 1		
Defined immunity to interference: Suitable for use in industrial areas		
Interference emission: Class		B

6 Accessories and Consumables

6.1 Accessories

This table contains an excerpt of the accessories that can be ordered. For information on other products, contact Sartorius.

	Product	Quantity	Order number
	Filter adapter for gelatine membrane filters	1	17801
	Tri-Clamp (TC 50.5) to bayonet adapter, stainless steel	1	17659---003

	Product	Quantity	Order number
	Bayonet to Tri-Clamp (TC 50.5) adapter	1	1ZAD-0025
	TC 50.5 to single-level hose nozzle DN 25 adapter, stainless steel. Can be used at the air outlet to drain the exhaust air.	1	17016
	Clamp for 1" to 1 1/2" sanitary flange	1	17033

6.2 Consumables

This table contains an excerpt of the consumables that can be ordered.
For information on other products, contact Sartorius.

	Product	Quantity	Order number
	Seal for Tri-Clamp connection TC 50.5, silicone	2	7EDSCV0003
	Sartofluor® Capsule air filter, 0.2 µm, 0.2 m² filter area	4	5181307T9--SS--A
	PVC hose, 2 m	1	17085
	PVC hose, 5 m	1	17088
	Silicone hose, 1 m	1	17662

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Phone: +49.551.308.0
Fax: +49.551.308.3289
www.sartorius.com

The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Masculine or feminine forms are used to facilitate legibility in these instructions and always simultaneously denote the other gender as well.

Copyright notice:

This instruction manual, including all of its components, is protected by copyright.

Any use beyond the limits of the copyright law is not permitted without our approval.

This applies in particular to reprinting, translation and editing irrespective of the type of media used.

© Sartorius Germany

Last updated:

03 | 2017

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Sartorius Lab Instruments GmbH & Co. KG